

Proline® Plus Mechanical Pipettes Dependable Durability



turning science **into solutions**

■ Proline® Plus Mechanical Pipettes

Dependable Durability

Sartorius' mechanical pipette Proline Plus in the Biohit *family* is designed to offer comfort and quality in every day manual pipetting. While it shares many of the excellent features of mLINE, it has a personal design and a robust feel for heavier use. In addition, it has the widest pipette range, including fixed volume pipettes, for when volumes need to be ready-set to avoid errors. Proline Plus pipettes are an excellent choice for both experienced laboratory professionals and students.

- The widest range of pipettes for various users and applications
- Ergonomic design with light pipetting forces, comfortable handle and finger support for reduced risk of strain injuries
- Minimised risk of contamination with Safe-Cone Filters and full autoclavability
- Highly durable with strengthened structure also for heavier use

Widest range of pipettes, both adjustable and fixed volume, available

- Adjustable volume pipettes for professionals
- Fixed volume pipettes for specific applications to avoid volume errors



Optiload tip loading mechanism in multichannel models

- Allows tip loading with an equal constant force onto every channel
- Enables perfect tip sealing onto every individual tip cone
- Makes tip ejection light and easy



□ Features and Benefits

Volume setting with click stop
– Easy to use volume setting

Ergonomic finger support
– No need to squeeze the pipette during pipetting, reduces the risk of RSI



Light pipetting forces & ergonomic design
– Prevents repetitive strain injury (RSI)

Volume colour-coding
– Eases correct tip volume selection through matching colour code

Thermal insulation of internal components
– Improves accuracy

Large numbers in volume display
– Quick and easy-to-read volume minimises risk of volume setting errors

Easy to use "do-it-yourself" calibration
– Quick re-calibration possible e.g. for different liquids

Materials have high chemical and UV-resistance
– Highly durable and long lasting pipette

Only three parts to disassemble
– The easiest pipette to clean and maintain, improves purity

Fully autoclavable without disassembling
– Easy to steam sterilise, improves purity



Safe-Cone Filters available for models > 10 µl
– Contamination free pipetting



**Widest Range of Pipettes**

The Proline Plus family includes the widest range of pipettes, from 3 µl | 5µl to 10 ml, as both adjustable and fixed volume devices. The fixed volume pipettes are particularly useful in student laboratories or in any laboratory where the application requires a constant volume, ready-set, to avoid errors.

Ergonomic Design

A comfortable handle and ergonomic finger support enable effortless pipetting. There is no need to squeeze the handle during pipetting, thereby reducing the risk of RSI.

**Perfect Tip Sealing with Optiload Tip Loading Mechanism in Multichannel Pipettes**

All Proline Plus multichannel pipettes are equipped with spring loaded tip cones – the Optiload mechanism secures even tip sealing and simultaneously allows tips to be loaded and ejected with minimum force. This is a particular advantage in multichannel models, where tip loading and ejecting generally require more force. Using less force reduces the risk of hand injuries.

**Prevent Contamination with Safe-Cone Filters**

Replaceable Safe-Cone Filters act as a barrier, preventing aerosols or fluids from reaching the internal components of the pipette. Safe-Cone Filters are available for all Proline Plus models greater than 10 µl. These filters must be replaced regularly, and in every case of over-aspiration.

Maintenance and Calibration Made Easy

No opening tools are needed when cleaning and maintaining Proline Plus pipettes, and only three parts need to be cleaned. These pipettes can be easily calibrated using the tool provided with them.

 Ordering Information

Proline® Plus

Cat. No.	Channels	Volume Range (µl)	Increment (µl)	Test Volume (µl)	Inacc. (%)	Impr. (%)	Safe-Cone Standard	Filters Plus	Optifit Tips (µl)	SafetySpace™ Filter Tips (µl)
728010	1	● 0.1–3	0.002	3	1.30	0.80	-	-	10	10
				1.5	2.40	1.60				
				0.3	10.00	6.00				
728020	1	● 0.5–10	0.01	10	1.00	0.60	-	-	10	10
				5	1.50	1.00				
				1	2.50	1.50				
728030	1	● 2–20	0.02	20	0.90	0.40	721014	-	200	20
				10	1.20	1.00				
				2	3.00	2.00				
728040	1	● 5–50	0.10	50	1.00	0.30	721008	721018	200	120
				25	1.40	0.50				
				5	3.00	1.50				
728050	1	● 10–100	0.10	100	0.80	0.15	721008	721018	200	120
				50	1.00	0.40			350	
				10	2.00	1.00				
728060	1	● 20–200	0.20	200	0.60	0.15	721007	721017	200	200
				100	0.80	0.30			350	300
				20	2.00	0.80				
728070	1	● 100–1000	1.00	1000	0.60	0.20	721006	721016	1000	1000
				500	0.60	0.20				
				100	1.00	0.40				
728080	1	● 500–5000	10.0	5000	0.50	0.20	721005	721015	5000	-
				2500	0.60	0.30				
				500	2.00	0.60				
728090	1	● 1–10 ml	20.0	10000	0.60	0.20	721005	721015	10 000	-
				5000	1.20	0.30				
				1000	3.00	0.60				
728120	8	● 0.5–10	0.01	10	1.50	1.00	-	-	10	10
				5	2.50	2.50				
				1	4.00	4.00				
728130	8	● 10–100	0.10	100	0.70	0.25	721008	721018	200	120
				50	1.00	0.70			350	
				10	3.00	1.50				
728140	8	● 30–300	0.20	300	0.60	0.25	721007	721017	350	300
				150	1.00	0.50				
				30	2.00	1.00				
728220	12	● 0.5–10	0.01	10	1.50	1.00	-	-	10	10
				5	2.50	2.50				
				1	4.00	4.00				
728230	12	● 10–100	0.10	100	0.70	0.25	721008	721018	200	120
				50	1.00	0.70			350	
				10	3.00	1.50				
728240	12	● 30–300	0.20	300	0.60	0.25	721007	721017	350	300
				150	1.00	0.50				
				30	2.00	1.00				

Ordering Information

Proline® Plus FIXED Volume, Single-Channel

Cat. No.	Channels	Volume (µl)	Increment (µl)	Test Volume (µl)	Inacc. (%)	Impr. (%)	Safe-Cone Standard	Filters Plus	Optifit Tips (µl)	SafetySpace™ Filter Tips (µl)
728515	1	5	–	5	1.30	1.20	–	–	10	10
728520	1	10	–	10	0.80	0.80	–	–	10	10
728530	1	20	–	20	0.60	0.50	721014	–	200	20
728535	1	25	–	25	0.50	0.30	721008	721018	200	120
728545	1	50	–	50	0.50	0.30	721008	721018	200	120
728550	1	100	–	100	0.50	0.30	721008	721018	200 350	120
728560	1	200	–	200	0.40	0.20	721007	721017	200 350	200 300
728565	1	250	–	250	0.40	0.20	721006	721016	1000	500 1000
728567	1	500	–	500	0.30	0.20	721006	721016	1000	500 1000
728570	1	1000	–	1000	0.30	0.20	721006	721016	1000	1000
728575	1	2000	–	2000	0.30	0.15	721005	721015	5000	–
728580	1	5000	–	5000	0.30	0.15	721005	721015	5000	–
728590	1	10 ml	–	10000	0.60	0.20	721005	721015	10 000	–



Carousel Stand

Ordering Information

Pipette Stands and Accessories

Cat. No.	Item
725600	Carousel Stand for 6 pipettes
725610	Pipette Holder for one pipette
725620	Linear Stand
726203	Calibration Tool Tube opener



Pipette Holder



Linear Stand



Calibration tool



□ Proline® Plus Starter Kits

Proline Plus Starter Kit offers an opportunity to test and get started with Proline Plus. You can choose between four Starter Kits, with two or three single-channel Proline Plus pipettes.

All kits include a range of useful accessories, such as a pipette holder and a calibration tool that also acts as a tube opener, as well as the required instructions.



□ Ordering Information

Cat. No.	Item
728650	Proline® Plus Starter Kit 1 Proline Plus 0.1 – 3 µl Proline Plus 0.5 – 10 µl Optifit Tip 0.1 – 10 µl, Single Tray × 2 Pipette Holder × 2 Calibration Tool Tube Opener × 2 Literature
728651	Proline® Plus Starter Kit 2 Proline Plus 0.5 – 10 µl Proline Plus 10 – 100 µl Proline Plus 100 – 1000 µl Optifit Tip 0.1 – 10 µl, Single Tray Optifit Tip 0.5 – 200 µl, Single Tray Optifit Tip 10 – 1000 µl, Single Tray Pipette Holder × 3 Calibration Tool Tube Opener × 3 Literature
728652	Proline® Plus Starter Kit 3 Proline Plus 2 – 20 µl Proline Plus 20 – 200 µl Proline Plus 100 – 1000 µl Optifit Tip 0.5 – 200 µl, Single Tray × 2 Optifit Tip 10 – 1000 µl, Single Tray Pipette Holder × 3 Calibration Tool Tube Opener × 3 Literature
728653	Proline® Plus Starter Kit 4 Proline Plus 500 – 5000 µl Proline Plus 1000 – 10000 µl Optifit Tip 100 – 5000 µl, Bulk Optifit Tip 1 – 10 ml, Bulk Pipette Holder × 2 Calibration Tool Tube Opener × 2 Literature



Contact Details

Sartorius Malaysia Sdn Bhd
L3-E-3B Enterprise 4
Technology Park Malaysia
Bukit Jalil 57000
Kuala Lumpur Malaysia

Phone +603.8996.0622

Fax +603.8996.0755

enquiry.malaysia@sartorius.com.my

Sartorius Biohit Liquid Handling Oy
Laippatie 1
FI-00880 Helsinki
Finland

Phone +358.9.755.951

lhinfo.finland@sartorius.com

www.sartorius.com

The status of the information, specifications and illustrations in this catalogue is indicated by the date given below.

Sartorius reserves the right to make changes to the technology, features, specifications and design of the equipment without notice.

All trademarks are Sartorius property unless otherwise stated. Patents granted or pending.



◀ www.sartorius.com