


Low EEO Agarose, Multi-purpose, Powder & Tablets



Gene Xpress LE Grade Agarose is a low EEO, multi-purpose, standard melting point agarose that yields high resolution sharp DNA bands with high clarity and low background. Its optimized gel strength enhances ease of gel processing and handling. It is manufactured using innovative Organic Solvent Free Manufacturing process that is greener and more environment friendly.

Application

Ideal for routine DNA and RNA gel electrophoresis and blotting.



0.75%, 1.0% and 1.25% Gene Xpress LE Agarose Gels,
1X TAE buffer
Lane A: EZ 200bp DNA Ladder (0.5ug/lane)
Lane B: EZ 1kb Plus DNA Ladder (0.5ug/lane)
5 V/cm for 2 hours

New!!
Only at
RM785.00*
per 500g



Instant LE Agarose Tablets Description

Instant Agarose Tablet, low EEO (LE), multi-purpose agarose in tablet form and delivered in blister convenient pack. Made of standard melting point agarose yields high resolution sharp DNA bands with high clarity and low background. Its optimized gel strength enhances ease of gel processing and handling.

Brand	Cat.No.	Product	Size
Gene Xpress	GN9051-500	LE Grade Agarose, Multi Purpose, Powder	500 g / btl
	GLE9012TAB-200	LE Grade Agarose, Multi Purpose, 0.5g Tablets	200 tablets / bx

Gene Xpress PLT (LLP0001852-LGN)

Contact Person: Chelston Cheah

Email: info@gene-xpress.com.my

Tel/Fax: +603-5740 8161 H/P: +6012-466 6885

