

AN EXCLUSIVE FACESHIELD



Anti-splitting splash

360° enclosed, effectively prevents from saliva splashing and fully protects the face and eyes.




Anti-dust

Effectively blocks sand, dust and fine particle.




Anti-fogging


Anti-fog coating on both sides of the shield can effectively prevent fog and maintain good visibility.



Front View



Side View



TEX-SHIELD SPECIFICATION

| | | |
|---|--------------------------|-------|
| 1 | APET Clear PlasticShield | 1 pcs |
| 2 | Foam | 1 pcs |
| 3 | Elastic Band | 1 pcs |
| 4 | Bonding Point | 2 |
| 5 | Shield Label | 1 pcs |

TEX-SHIELD PRODUCT DESCRIPTION

Non-Sterile, Single Use Face Shield as a physical barrier for face from splashing and spraying. Clear Face Shield Film, Elastic Band and Comfortable Foam.

Size : 330mm x 250mm
Material : Non-Sterile +Foam +Elastic Band
Usage Range : Suitable for use in working environments that require face protection.

Note:

1. When wearing Tex-Shield, ensure it is securely fasten around the forehead. Avoid touching on plastic sheet to minimize stain that can cause visual obstruction.
2. Tex-Shield shall prevents direct contact against non-hazardous liquid direct splash or spraying. Recommended to be worn with face mask for better protection.
3. Tex-Shield is designed as a single use product, discard if damage. It is not design to used for fumes, mist, gases, vapor and impact protection.
4. Tex-Shield is recommended be discarded as medical waste when contaminated with blood or bodily fluid.

Distributed by:



Biohaus Asia Sdn Bhd
 No. 2A, Jalan Sepadu 25/123A,
 Axis Industrial Park, Seksyen 25,
 40400 Shah Alam
 estore.biohaus.my

