FAVORIT® EQUIPMENTOVERVIEW

WATER STILL

FAVORIT W4L

Cat. No.: P21001



FEATURES

- Horizontal thick wall heat resistant borosilicate glass boiler fitted with 3kW chromium plated heater.
- High efficiency coil condenser ensure low temperature, high output quality distillate of conductivity 3 to 4µs/cm.
- Compact construction suitable for bench and wall mounting.
- Screwthread connectors enable easy fitting or removal of hose safely.
- Built-in constant level device with Vendura-Nova valve for drainage provides easy descaling without dismantling.
- Safety thermostat built into the heater protect the Favorit water still in the event of water supply failure and automatically resets when the water supply is restored.

SPECIFICATIONS

Output: litres/hour	4L/hour
Quality of distillate:	
ρН	5.0 - 6.5
Conductivity µS/cm	3.0 - 4.0
Resistivity megOhm-cm	0.25 - 0.30
Pyrogen content	Pyrogen free
Temperature	35°C-50°C

Water Supply	1 litre / min Min. pressure 3 psi (0.2 x 10 ⁵ NM ⁻²)
Power	220 / 240V 50 / 60 Hz single phase 13A fused
Dimensions : (w x d x h)	50 x 15 x 45cm
Catalogue Number:	P21001