



Key Features

- Reliable Performance
- Energy Saving
- Low Noise

More Detail:

http://www.haiermedical.com/products/preservation/86_ult_freezer/dw-86l486.html

Specifications

Technical Data	
Type	-86°C ULT Freezer
Cabinet Type	Upright
Climate Class	N
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	HFC
Noise Level (dB)	49
Performance	
Cooling performance(°C)	-86
Temperature Range(°C)	-40~-86
Control	
Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50/208-230/60/115/60		
Power(W)	1000	1000	1200
Electrical Current(A)	9	5.5	18

Dimensions

Capacity(L/Cu.Ft)	486/17.1		
Net/Gross Weight(approx)	290/310(kg) 639.3/683.4(lbs)		
Interior Dimensions(W*D*H)	590*630*1310(mm) 23.3*24.8*51.6(in)		
Exterior Dimensions(W*D*H)	945*900*1980(mm) 37.2*35.4*78.0(in)		
Packing Dimensions(W*D*H)	995*995*2150(mm) 39.2*39.2*84.6(in)		
Container load (20'/40'/40'H)	12/24/24		

Functions

High/Low Temperature	Y
Remote Alarm	Y
Hot Condenser	Y
Power Failure	Y
High/Low Voltage	Y
Sensor Error	Y
Low Battery	Y
High Ambient Temp	Y
Door Ajar	Y

Accessories

Caster	Y
Foot	Y
Test Hole	Y
Shelves/Inner Doors	3/4
Temperature Recorder	Optional
USB Interface	Y

Rs232/485 Interface	Optional
CO2 Backup System	Optional
LN2 Backup System	Optional
Racks And Boxes	Optional

CO2 Backup System

LN2 Backup System

Chart Recorder

Others

Certificate	CE/UL
--------------------	-------