



Key Features

- Reliable Performance
- Energy Saving
- Low Noise

More Detail:

http://www.haiermedical.com/products/preservation/86_ult_freezer/dw-86l628.html

Specifications

Technical Data

Type -86°C ULT Freeer

Cabinet Type Upright

Climate Class N

Cooling Type Direct cooling

Defrost Mode Manual

Refrigerant HFC

Noise Lever (dB) 49

Performance

Cooling performance(°C) -86

Temperature Range(°C) -40~-86

Control

Controller Microprocessor

Display LED

Electrical Data

Power Supply(V/Hz) 220-240/50/208-230/60

Power(W)	1200	1200
Electrical Current(A)	7	9
Dimensions		
Capacity(L/Cu.Ft)	628/22.2	
Net/Gross Weight(approx)	301/323(kg) 664.0/712.0(lbs)	
Interior Dimensions(W*D*H)	760*630*1310(mm) 29.9*24.8*51.6(in)	
Exterior Dimensions(W*D*H)	1035*900*1980(mm) 40.7*35.4*78.0(in)	
Packing Dimensions(W*D*H)	1080*965*2150(mm) 42.5*38.0*84.6(in)	
Container load (20'/40'/40'H)	12/24/24	
Functions		
High/Low Temperature	Y	
Remote Alarm	Y	
Hot Condenser	Y	
Power Failure	Y	
High/Low Voltage	Y	
Sensor Error	Y	
Low Battery	Y	
High Ambient Temp	Y	
Door Ajar	Y	
Accessories		
Caster	Y	
Foot	Y	
Test Hole	Y	
Shelves/Inner Doors	3/4	
Temperature Recorder	Optional	
USB Interface	Y	
Rs232/485 Interface	Optional	
CO2 Backup System	Optional	
LN2 Backup System	Optional	
Racks And Boxes	Optional	

Racks & Boxes(Optional)

Rack Model	Rack Description	Rack Dimensions (W*D*H)	Q'ty of Racks	Q'ty of Boxes	Q'ty of Samples 9*9 box / 10*10 box
DCJ-54-A	2 inch stainless steel rack 5*4	139*544*287(mm) 5.5*21.4*11.3(in)	20	400	32400/40000
DCJ-54-B	2 inch stainless steel drawer rack 5*4	140*554*289(mm) 5.5*21.8*11.4(in)	20	400	32400/40000

CO2 Backup System

LN2 Backup System

Chart Recorder

Others

Certificate

CE/UL