

# Pharmacy Refrigerator

## — HYC-940



### Six Alarm Functions

High/low temperature, remote alarm, power failure, sensor error, low battery, door ajar



### Air Cooling, Auto Defrost

Better cabinet temperature uniformity, auto defrost



### Adjustable Shelves

The shelves are adjustable to any height to fulfill different requirements



### Heated Double Pane Glass Door

Clear display of stored products without condensation at 32°C and 85% relative humidity condition



### Self Closing Door

Energy conserving, maintains temperature uniformity to protect valuable samples

## Ergonomic Design



High Accuracy  
Temperature Control



Interior  
LED



Lockable  
Design

## Qingdao Haier Biomedical Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone,  
Qingdao, 266109, P.R. China  
Tel: +86-0532-88935593  
Website: [www.haiermedical.com](http://www.haiermedical.com)



Haier Biomedical  
International



Haier Biomedical  
International



@haiermedicalint



@haierbiomedicalint



Haier Biomedical  
International

## Specifications



### Technical Data

Cabinet Type	Upright, double glass door
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	HC
Noise((dB(A))	43

### Performance

Temp Range(°C)	2~8
----------------	-----

### Control

Controller	Microprocessor
Display	LED

### Electrical Data

Power Supply(V/Hz)	220~240/50/60	115/60
Power(W)	700	850
Electrical Current(A)	3.5	7.5

### Dimensions

Capacity(L/Cu.Ft)	890/31.4
Net/Gross Weight(approx)	227/279(kg) 500.4/615.1(lbs)
Interior Dimensions(W*D*H)	1030*590*1425(mm) 40.6*23.2*56.1(in)
Exterior Dimensions(W*D*H)	1130*755*1980(mm) 44.5*29.7*78.0(in)
Packing Dimensions(W*D*H)	1250*850*2130(mm) 49.2*33.5*83.9(in)
Container Load(20'/40'/40'H)	8/18/18



### Functions

High/Low Temp	Y
Remote Alarm	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

### Accessories

Caster	Y
Foot	Y
Porthole	Y
Shelves/Basket	12/-
Temperature Recorder	Optional
USB Interface	Optional

### Others

Certificate	CE	UL
-------------	----	----