



# Cotton Tipped Applicator

Wood handle



## Description

This absorbent cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a rigid white birch handle by an aqueous based adhesive.

## Dimensions

**6" Overall Length (152.4 mm)**

### Tip

- 0.188" Dia. (4.775 mm)
- 0.625" Length (15.875 mm)

### Handle

- 0.1" Dia. (2.54 mm)
- 5.812" Length (147.624 mm)

## Qualities

- Medical grade quality
- Standard size cotton tip
- Absorbent tip
- Soft smooth tip
- Made in USA

## Applications

- Specimen collection
- Wound care
- Cell collection for DNA testing
- Medication applicator
- First-aid kit component
- Site preparation for health/beauty and cosmetic applications
- Tattoos

Ref #	Packaging
25-806 1 WC Sterile	10/100/1 (10 boxes of 100 individually wrapped packages, 1 p/pkg)
25-806 2WC Sterile	10/100/2 (10 boxes of 100 individually wrapped packages, 2 p/pkg)
25-806 5WC Sterile	10/20/5 (10 boxes of 20 individually wrapped packages, 5 p/pkg)
25-806 10WC Sterile	10/20/10 (10 boxes of 20 individually wrapped packages, 10 p/pkg)



**Puritan®**