

# Oil Free Vacuum Pump



## Features

← **No air pollution, maintenance free**

Rocker series pumps are driven by piston, without the need of lubricant, regular oil changes and maintenance; with no oil pollution.

← **Prote™ overflow protection**

Rocker series pumps are equipped with Prote™ overflow protection which can prevent a large amount of liquid from being sucked into the pump.

← **Moisture trap with filter cartridge**

Rocker series pumps are equipped with filter cartridge at the air inlet to filter particle and moisture to prolong the life of pump.

← **Vacuum regulator**

Rocker series pumps are equipped with vacuum regulator to adjust vacuum.

← **Thermal protection device**

Every motor of Rocker series pumps has a built-in thermal protection device to shut off the pump automatically when overheated and then resume working when the temperature cools down.

← **International certification**

CE certification  
 CSA certification (Rocker 300)  
 IP 30 protection (Rocker 300)  
 RoHS certification

← **Two-Year Warranty**

## Applications

- ← Vacuum filtration
- ← Suspended solid test
- ← Microbiological test
- ← Air sampling

## Ordering Information

**167300-11(22)**

Rocker 300, Vacuum Pump  
 AC110V, 60Hz (AC220V, 50Hz)

**167330-00**

Rocker 300DC, DC Power Vacuum Pump (DC12V)



## Specification

Model	Rocker 300		Rocker 300DC
Voltage / Freq.	110V/60Hz	220V/50Hz	DC12V
Max. power	65W	60W	35W
Max. current	0.7A	0.3A	3A
Max. vacuum	105 mbar		145 mbar
Max. flow rate	23 L/min	20 L/min	25 L/min
Motor rotation speed	1750 RPM	1450 RPM	1800 RPM
Horse power	1/8 HP		1/8 HP
Noise level	50 dB		50 dB
Hose barb	ID8		ID8
Vacuum regulator	Yes		Yes
Overflow protection/ Filter cartridge	Yes		Yes
Safety	Thermal Protection		Thermal Protection
Net weight	4.1 Kg		3.5 Kg
Dimension (LxWxH)	26.8 x 13.5 x 20.4 cm		27.2 x 12.1 x 16.5 cm

# Pump Accessories

## Hose Barb



Cat. No.:	167300-20	167300-21	167300-63	167300-23	167300-24	187100-32	187100-31	187100-34	167300-25
ID:	3/8 inch	5/16 inch	1/4 inch	3/16 inch	1/8 inch	1/2 inch	3/8 inch	1/4 inch	5/16 inch
	ID10	ID8	ID6	ID5	ID3	ID12	ID10	ID6	ID8
Thread:	1/8 PT	1/8 PT	1/8 PT	1/8 PT	1/8 PT	1/4 SAE	1/4 SAE	1/4 SAE	1/8 PT

## Moisture Trap / Vacuum Regulator / Pressure Regulator



167301-06	167801-06	169411-06	169611-06	189300-06	167320-06	167440-06
For Rocker 300/400/ 410/430	For Rocker 801/811	For Chemker 410/411	For Chemker 610/611	For Rocker 300C/400C/410C	For Rocker 320/420	For Rocker 440

## Filter Cartridge



167301-07	167801-07
For Rocker 300/400/410/430 (with Prote™ overflow protection)	For Rocker 801/811 (without Prote™ overflow protection)

## Chemical Resistant Hose Barb



169300-21	169300-20
5/16 inch ID8 1/4 NPS	3/8 inch ID10 1/4 NPS

## Air Outlet (Copper Muffler)



167300-42	167600-42
For Rocker 300/400/410/430	For Rocker 800/801/810/811

## Check Valve



187100-33
For Tanker 130/150/215/230 ID10 1/4 SAE

## Silicone Tube



167200-38	199100-68	180300-68
Ø 6 x 12 mm Length: 1 m	Ø 4.2 x 8 mm Length: 2 m	Ø 8 x 14 mm Length: 2 m

## Universal Foot Switch



167200-41	167200-41-1
Momentary action	Maintained action